

9-POINT OF PAIN PROTOCOL

Avazzia Pro Sport

Point of Pain
(Counterclockwise)

Contra Lateral
(Clockwise)

IR □ 3	IR □ 2	IR □ 9
IR □ 4	IR/D □ 1	IR □ 8
IR □ 5	IR □ 6	IR □ 7

IR □ 9	IR □ 2	IR □ 3
IR □ 8	IR/D □ 1	IR □ 4
IR □ 7	IR □ 6	IR □ 5

- Turn on your device and power up to a comfortable pricking sensation leaving the device on default (Relax Asses). Press the left arrow to switch to Biofeedback mode.
*Note - Remember to record the readings of your IR (Initial Reaction, top center display), D (Dose, ongoing reaction bottom left corner) and Z (Zero, ongoing reaction bottom left corner) for comparison.
- Locate Point of Pain (POP). This will be your first spot; you will always bring this location to a D (dose).
*Note- Before bringing a spot to D, always remove device from skin and replace on same location.
- Proceed with IRs in the order of the numbered boxes (counterclockwise) and only D the locations that have an IR that is higher than your initial Point of Pain IR.
- Go to the contra lateral side of point of pain and do the same procedure, only clockwise.
- After you have IR/D your POP and Contra lateral sides, you will find the single highest D from both sides and take that spot to Z.
- Next, we incorporate the CNS. Remove device from the body and place it on the position of the Zero. Move the device horizontally until you reach the corresponding position at the center of the spine. Take an IR reading, remove and replace to Dose this spot.